

Bänkvåg LEM7
Bench scale LEM7

Manual/Förenklad handhavande

LEM7



DAGLIG ANVÄNDNING:

Uppstart och vägning utan tara:

1. Placera vågen på stabilt underlag och i "våg".
2. Slå på huvudströmbrytaren. Finns på höger sida under vågen.
3. Starta sedan vågen via "ON/OFF" tangent.
4. Är det information i viktdisplayen tryck "ZERO" tangent.
5. Påbörja vägning.

Vägning med tara:

6. Placera behållare på vågplattan och tryck "ZERO" tangent.
7. Påbörja vägning.
8. Ta bort tara genom att trycka på "ZERO" när vågplattan är tom.

Växla mellan viktenheter:

9. Tryck på " UNITS " för att byta vikt enhet.
(Det finns 14 st. viktenheter man kan välja mellan.)

Avstängning av batteri-spar funktionen:

10. Håll " ZERO " intryckt när vågen startas.
A-ON visas i displayen.
11. Tryck på " ZERO " en gång.
A-OFF visas i displayen.
12. Stäng av vågen och slå på den för att spara och återgå till vägning.

Styckräkning:

13. Lägg på önskat antal för referens (10-20-50 eller 100st).
14. Tryck på " PCS " (S-10pcs visas i displayen).
15. Tryck på " UNITS " för att välja rätt antal.
16. Tryck sedan på " PCS " för att bekräfta.
17. För att återgå till vägningssläge tryck på " UNITS ".

DAILY USE:

Start and weighing without tare:

1. Place the scale on a stable bench or other surface.
2. Switch on the main power switch at the bottom right of the scale.
3. Press **“ON/OFF”**
4. Press **“ZERO”** if needed to zero the scale
5. Start weighing

Weighing with tare

1. Place the container on the scale platform and press **“ZERO”**
2. Start weighing
3. End weighing with tara by pressing **“ZERO”** with an empty platform.

Change between units

- Press **“UNITS”** to change between different weighing units (There are 14 units to choose between)

Turn off battery save mode

1. Press **“ZERO”** when starting the scale
A-ON is displayed
2. Press **“ZERO”** again
A-ON is displayed
3. Tur off the scale and start it again. The scale is ready to use.

Counting

1. Place a sample on the scale platform (10-20-50 or 100 units).
2. Press **“PCS”**. (S-10pcs is displayed).
3. Press **“UNITS”** to decide the correct number of samples.
4. Press **“PCS”** to confirm.
5. To exit counting mode: press **“UNITS”**.